



HELEN

HELEN HOLDING GROUP

HELEN

HELEN HOLDING GROUP



Departure Destination | Saudi Arabia

FROM 2024



Destination Reaching | The World

TO 2025

SMART SKIN ANA

WWW.HELENHOLDING.COM - 920033794 - HELENHOLDING@OUTLOOK.SA



9 2 0 0 3 3 7 9 4

HELEN SKIN ANALYSIS

Magic Mirror Max

Magic Mirror Max-Skin Analysis System is the world's most advanced skin detection equipment together with shooting, analyzing, displaying 3 in 1. It adopts RGB, UV, PL spectral imaging technology, combines with artificial intelligence and image analysis, 12 years market testing, 30 million clinical database, achieving 15 seconds efficient skin analysis. Such as: spots, pores, moisture, textures, wrinkles and so on 12 epidermis and epidermis skin problems. According to the test results, Magic Mirror Max-recommend products for tailor-made skin management program. It is an essential skin detection equipment for cosmetic company, beauty salons, skin clinics, training schools, the market demand rate as high as 100%.

for



easy to
perform
very
effective





WORKING WAY !!

Based on big data, combined with the latest Deep Learning technology, provides accurate facial analysis area, automatically identify the analysis area.

Combining artificial intelligence and the most sophisticated graphics analysis algorithms, makes the analysis more scientific.

Helen


MEDICAL HOLDING COMPAN

Helen



Helen

MEDICAL HOLDING COMPAN



Upoo

RGB Pore: represents current big pores on surface skin which is bigger than normal size pore over diameter 0.02-0.05 mm.

RGB Spot: represents current spot which color is darker than normal skin color and with round or irregular round shape on surface skin.


RGB Wrinkle : represents current wrinkle condition as well as area which is not smooth and not flat, the higher the percentage, the better of wrinkle condition.




Helen

MEDICAL HOLDING COMPAN


OVER VIEW

VISION





1 general test +12 single test report+SD card+PDF print



HELEN SKIN ANALYSIS

Skin analysis

Specifications

Power	45W	Machine Size	42 * 36 * 56cm
Voltage	110 ~ 230 VAC $\pm 10\%$	Package Size	50 * 45 * 67cm
Current	0.2A 50HZ	Package Method	Carton Box
Camera Pixel	20 Megapixel	Net Weight	8 KG
Screen Size	10.1 inch	Gross Weight	12 KG



WWW.HELENHOLDING.COM

HELENHOLDING@OUTLOOK.SA

KAS- ABHA – ALMAHALA DIS - TEL | 920033794



HELEN

HELEN HOLDING GROUP